

P/N 961-01479-07 Article Structure Description	
Liner Type	Split
Size	ID 3.5mm; OD 5mm; length 95mm
Injector	SPL-2010 Plus, SPL-2010
Packing	Without packing
Volume (µl)	860
Deactivation	Siltek-deactivated
Max Temperature	450°C
Remarks	Split liner without packing material (no glass wool or CarboFrit). Siltek-deactivated for excellent inertness for a wide sample range, difficult matrixes and reactive compounds. Ideal for chlorinated pesticide analysis, not recommended for basic compounds.
Packing	5 liner per package

