

<b>P/N 961-01480-02 Article Structure Description</b>	
Liner Type	Splitless
Size	ID 3.5mm; OD 5mm; length 95mm
Injector	SPL-2010 Plus, SPL-2010, SPL-2014 WBI-2010 Plus, WBI-2010, WBI-2014, SPL-2025
Packing	Glass Wool
Volume (µl)	670
Deactivation	Phenylmethyl-deactivated
Max Temperature	350°C
Remarks	Phenylmethyl-deactivated (IP) for better recovery of polar and nonpolar compounds.
Packing	5 liners per package

