

P/N 961-01480-05 Article Structure Description	
Liner Type	Splitless
Size	ID 3.5mm; OD 5mm; length 95mm
Injector	SPL-2010 Plus, SPL-2010, SPL-2014 WBI-2010 Plus, WBI-2010, WBI-2014, SPL-2025
Packing	Glass Wool
Volume (µl)	670
Deactivation	Base-deactivated
Max Temperature	350°C
Remarks	Base-deactivated for better recovery of basic compounds like amines.
Packing	5 liner per package

