

| P/N 961-01480-06 Article Structure Description | |
|--|---|
| Liner Type | Splitless |
| Size | ID 3.5mm; OD 5mm; length 95mm |
| Injector | SPL-2010 Plus, SPL-2010, SPL-2014 WBI-2010 Plus, WBI-2010, WBI-2014, SPL-2025 |
| Packing | CarboFrit |
| Volume (µl) | 670 |
| Deactivation | Base-deactivated |
| Max Temperature | 350°C |
| Remarks | Base-deactivated for better recovery of basic compounds like amines. |
| Packing | 5 liner per package |

