

P/N 961-01480-08 Article Structure Description	
Liner Type	Splitless
Size	ID 3.5mm; OD 5mm; length 95mm
Injector	SPL-2010 Plus, SPL-2010, SPL-2014 WBI-2010 Plus, WBI-2010, WBI-2014, SPL-2025
Packing	Glass Wool
Volume (µl)	670
Deactivation	Siltek-deactivated
Max Temperature	450°C
Remarks	Siltek-deactivated for excellent inertness for a wide sample range, difficult matrixes and reactive compounds. Ideal for chlorinated pesticide analysis, not recommended for basic compounds.
Packing	5 liner per package

