

P/N 961-01481-08 Article Structure Description	
Liner Type	Split liner for PTV injector
Size	ID 1.5mm; OD 3.5mm; length 95mm
Injector	PTV-2010
Packing	Glass Wool
Volume (µl)	170µl
Deactivation	Siltek-deactivated
Max Temperature	450°C
Remarks	Siltek-deactivated for excellent inertness for a wide sample range, difficult matrixes and reactive compounds. Ideal for chlorinated pesticide analysis, not recommended for basic compounds.
Packing	5 liner per package

