

<b>P/N 961-01482-01 Article Structure Description</b>	
Liner Type	SPME liner for SPL injector
Size	ID 0.75mm; OD 5.0mm; length 95mm
Injector	SPL-2010, SPL-2010 Plus, SPL-2014, WBI-2010, WBI-2010Plus, WBI-2014, and SPL-2025
Packing	Without packing
Volume (µl)	42µl
Deactivation	Phenylmethyl-deactivated
Max Temperature	350°C
Remarks	Phenylmethyl-deactivated (IP) for better recovery of polar and nonpolar compounds. Solid Phase Micro Extraction (SPME) is a special sampling procedure. Best results are obtained with injection into liners with small volume for better sample transfer on the column. Liner packing materials are not needed with this technique.
Packing	5 liner per package

