

P/N 961-01482-07 Article Structure Description	
Liner Type	SPME liner for SPL injector
Size	ID 0.75mm; OD 5.0mm; length 95mm
Injector	SPL-2010, SPL-2010 Plus, SPL-2014, WBI-2010, WBI-2010Plus, WBI-2014, and SPL-2025
Packing	Without packing
Volume (µl)	42µl
Deactivation	Siltek-deactivated
Max Temperature	450°C
Remarks	<p>Siltek-deactivated for excellent inertness for a wide sample range, difficult matrixes and reactive compounds. Ideal for chlorinated pesticide analysis, not recommended for basic compounds.</p> <p>Solid Phase Micro Extraction (SPME) is a special sampling procedure. Best results are obtained with injection into liners with small volume for better sample transfer on the column. Liner packing materials are not needed with this technique.</p>
Packing	5 liner per package

